



TF Card



Air filter
x5



User manual



Nasal Mask



Air tubing



Power supply



Travel bag

Auto CPAP VM-6

American Sensor

German Technology

Intelligent Breathing Algorithm

Improve Overall Breathing

Whisper quiet , light and compact

Ultra-quiet design, noise <28 dB; Small size, light weight, and portable design.

Effective Treatment

AST(Automatic Synchronization Technology) anticipates respiratory waveform and provides corresponding pressure to ensure patient's comfort and mechanical synchrony.

Efficient heating humidifer system

Integrated thermostatic humidifier, Anti-backflow desgin, avoid dry burning or water overflow.

Professional sleep study report

Unique sleep study report of whole week & month,New operating system is easier to operate , interface is more exquisite.



Specification

Description	Auto CPAP VM-6
Work Mode	CPAP, Auto
Sound	<28 dB
Screen	2.8 inch TFT, 320 x 240
Pressure Range	4-20 cmH ₂ O, auto adjustment
Data Storage	TF Card
Power Supply	Input 100-240V, 50-60Hz, 1-2A; Output + 24V, 2.5A
DC Voltage	24VDC, 2.5A
Over-Temperature Protection	Stop automatically, alarm with 4 times
Lock Function	Yes
Screen Auto Sleep	Auto sleep in 15 minutes, press any button to highlight screen
Ramp Item	0-60min adjustable 1 min/Step
Auto ON/Off	Yes
EPR	Yes, 0-3 level
Auto-Altitude	Yes
Leak Compensation	Yes
Awake Detection	Yes
Breath Event Detection	Yes
Auto Pressure Titration	Yes
Water Capacity	200ML max
Humidifier Level	Integrated Heated, 0-5 level
Alert	Power off, mask off
Machine Dimensions(cm)	29L x 15W x 10H(machine)
Carton Dimensions(cm)	34.5L x 20.5W x 25.5H(1pc carton)
	40L x 34.5W x 25.5H(2pcs carton)
	59.5L x 34.5W x 25.5H(3pcs carton)
Weight(with humidifier)	1.7Kg(net weight); 4.0Kg(gross weight)
Work For	OSA



Hunan VentMed Medical Technology Co., Ltd.

Factory: Baoqing industrial park, Shaoyang, Hunan, China

Office: +86-739-5188959 Web: www.ventmed.net

